

© 2006 General Electric Company – All rights reserved.
GE and GE Monogram are trademarks of
General Electric Company.
GE Healthcare Oy., doing business as GE Healthcare.

Healthcare Re-imagined

GE is dedicated to helping you transform healthcare delivery by driving critical breakthroughs in biology and technology. Our expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, and biopharmaceutical manufacturing technologies is enabling healthcare professionals around the world discover new ways to predict, diagnose and treat disease earlier. We call this model of care “Early Health.” The goal: to help clinicians detect disease earlier, access more information and intervene earlier with more targeted treatments, so they can help their patients live their lives to the fullest. Re-think, Re-discover, Re-invent, Re-imagine.

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11 • Fax +358 9 146 3310
www.gehealthcare.com



GE Healthcare SpO₂ Connection Guide

GE Healthcare

SpO₂ Connection Guide



M1071180

	Sensor Compatibility	Monitor	Page No.				
2-13	Ohmeda OxyTip+	3775	4				
		37xx Series	4				
		3800	4				
14-21	Masimo LNOP	3900	4				
		3900P	4				
		4700 respiratory gas monitor	4				
22-27	Nellcor	5250 respiratory gas monitor	4				
		Capnomac Ultima	6	8			
		Cardiicap/5, Cardiicap	6	8			
		Cardiicap/5 with N-XOSAT	4				
30-31	Nellcor Compatibles	CardioLab®/Mac-Lab®					
		with Tram Module 451M	14	16	22	32	
		CardioServ	26	30			
32-33	Masimo LNCS	Corometrics® 118 & 556	26				
		Corometrics® 120	22				
		Corometrics® 2120	22				
		Corometrics® 259	10	14	16	22	32
		Dash 2000	22	30			
		Dash 3000	14	16	22	30	32
		Dash 4000	14	16	22	30	32
		Dash 5000	14	16	22	30	32
		DINAMAP™ Compact	26				
		DINAMAP™ MPS	26				
		DINAMAP™ Pro 1000 V3	18	20	24	26	32
		DINAMAP™ Procure	18	20	24	32	
		DINAMAP™ ProSeries V1	26	30			
		DINAMAP™ Pro 1000 V1 and V2	26	30			
		DINAMAP™ ProSeries V2	18	20	24	26	32
		E-NSAT Modules for i4	28				
		Eagle 1000/3000/4000 (rectangle connector)	22	30			
		Eagle 1000/3000/4000 (round connector)	26	30			
		AS/3, CS/3, S/5 Modular Monitors with M-series modules	6	8			
		with M-NSAT	28				
		with M-OSAT	8				
		S/5, i4 Modular Monitors with E-series modules	10				
		with E-NSAT	28				
		MRI Monitors	12				
		Ohmeda TruSat™	2				
		Ohmeda TuffSat™	4				
		Responder 3000	22				
		S/5 Light Monitor	6	8			
		Satlite Plus II	6	8			
		Solar 8000/8000M, i/9500 with Tram Module 451M/851M	14	16	22	30	32

Instructions For Use

For detailed information on the use of these products, always refer to the User Guide/Instructions For Use accompanying the product/device.

	Sensor Compatibility	Monitor	Page No.			
2-13	Ohmeda OxyTip+	3775	4			
		37xx Series	4			
		3800	4			
14-21	Masimo LNOP	3900	4			
		3900P	4			
		4700 respiratory gas monitor	4			
		5250 respiratory gas monitor	4			
		Capnomac Ultima	6	8		
22-27	Nellcor	Cardiocap/5, Cardiocap	6	8		
		Cardiocap/5 with N-XOSAT	4			
		CardioLab®/Mac-Lab® with Tram Module 451M	14	16	22	
		CardioServ	26			
		Corometrics® 118 & 556	26			
		Corometrics® 120	22			
		Corometrics® 2120	22			
		Corometrics® 259	10	14	16	22
		Dash 2000	22			
		Dash 3000	14	16	22	
		Dash 4000	14	16	22	
		Dash 5000	14	16	22	
		DINAMAP™ Compact	26			
		DINAMAP™ MPS	26			
		DINAMAP™ Pro 1000 V3	18	20	24	26
		DINAMAP™ Procure	18	20	24	
		DINAMAP™ ProSeries V1	26			
		DINAMAP™ Pro 1000 V1 and V2	26			
		DINAMAP™ ProSeries V2	18	20	24	26
		E-NSAT Modules for i4	28			
		Eagle 1000/3000/4000 (rectangle connector)	22			
		Eagle 1000/3000/4000 (round connector)	26			
AS/3, CS/3, S/5 Modular Monitors with M-series modules	6	8				
with M-NSAT	28					
with M-OSAT	8					
S/5, i4 Modular Monitors with E-series modules	10					
with E-NSAT	28					
MRI Monitors	12					
Ohmeda TruSat™	2					
Ohmeda TuffSat™	4					
Responder 3000	22					
S/5 Light Monitor	6	8				
Satlite Plus II	6	8				
Solar 8000/8000M, i/9500 with Tram Module 451M/851M	14	16	22			

Instructions For Use



For detailed information on the use of these products, always refer to the User Guide/Instructions For Use accompanying the product/device.

Monitor



Ohmeda TruSat™
Ohmeda TruSat™ TD

Connecting Cables



Sensor Type



OXY-F-UN Finger sensor with UN connector, 1 m, requires interconnect cable (1/box)



OXY-W-UN Wrap sensor with UN connector, 1 m, requires interconnect cable (1/box)



OXY-SE-3 Sensitive Skin sensor with UN connector, 4 m, requires interconnect cable (3/box)



OXY-AP-25 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (25/box)



OXY-AP-10 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (10/box)




OXY-AF-10 AllFit adhesive sensor, 0.9 m, extra adhesive tape included, for adult, pediatric and neonatal use, requires interconnect cable (10/box)



OXY-E-UN Ear sensor with UN connector, 1 m, requires interconnect cable (1/box)



OXY-F4-MC Finger sensor with MC connector, 4 m (1/box)



OXY-E4-MC Ear sensor with MC connector, 4 m, headband and clip included (1/box)

Accessories

- OXY-RTW** OxyTip+™ wide replacement tape, adhesive, for use with OxyTip+™ Wrap, OxyTip+™ Sensitive Skin sensor or OxyTip+™ AllFit 100/box
- OXY-RWL** Foam wrap replacement, large, weight range ≥3 kg, for use with OxyTip+™ Wrap sensor 24/box
- OXY-RWM** Foam Wrap replacement, medium, weight range ≥3 kg, for use with OxyTip+™ Sensitive Skin Sensor 24/box
- OXY-RWS** Foam Wrap replacement, small, weight range <3 kg or small finger, for use with OxyTip+™ Sensitive Skin sensor 24/box
- OXY-RTB** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, bears 100/box
- OXY-RT** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, blue 100/box
- OXY-SND** Infant foam sandal, for use with OxyTip+™ Sensitive Skin sensor 3/box

Monitors




Ohmeda TuffSat™
3775, 37xx Series, 3800, 3900,
3900P, 4700 and 5250 respiratory
gas monitors

Cardiocap/5 with N-XOSAT option
S/5, AS/3, CS/3 Modular Monitors
with M-OSAT module

Connecting Cables



Interconnect cable with "H" (Ohmeda)
connector and "UN" connector (1/box)
OXY-OL1 (1m)
OXY-OL3 (3m)



Integrated cable with "H"
connector and sensor

Sensor Type



OXY-F-UN Finger sensor with UN connector, 1 m,
requires interconnect cable (1/box)



OXY-W-UN Wrap sensor with UN connector, 1 m,
requires interconnect cable (1/box)



OXY-SE-3 Sensitive Skin sensor with UN connector,
4 m, requires interconnect cable (3/box)



OXY-AP-25 Adult/Pediatric adhesive sensor, for adult
and pediatric use, requires interconnect
cable (25/box)



OXY-AP-10 Adult/Pediatric adhesive sensor, for adult
and pediatric use, requires interconnect
cable (10/box)



OXY-AF-10 AllFit adhesive sensor, 0.9 m, extra adhesive
tape included, for adult, pediatric and neonatal
use, requires interconnect cable (10/box)



OXY-W4-H Wrap sensor with H (Ohmeda) connector,
1 m (1/box)



OXY-F4-H Finger sensor with H (Ohmeda) connector,
4 m (1/box)



OXY-F1-H Finger sensor with H (Ohmeda) connector,
1 m (1/box)



OXY-E4-H Ear sensor with H (Ohmeda) connector,
4 m (1/box)

Accessories

- OXY-RTW OxyTip+™ wide replacement tape, adhesive, for use with OxyTip+™ Wrap, OxyTip+™ Sensitive Skin sensor or OxyTip+™ AllFit 100/box
- OXY-RWL Foam wrap replacement, large, weight range ≥3 kg, for use with OxyTip+™ Wrap sensor 24/box
- OXY-RWM Foam Wrap replacement, medium, weight range ≥3 kg, for use with OxyTip+™ Sensitive Skin Sensor 24/box
- OXY-RWS Foam Wrap replacement, small, weight range <3 kg or small finger, for use with OxyTip+™ Sensitive Skin sensor 24/box
- OXY-RTB OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, bears 100/box
- OXY-RT OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, blue 100/box
- OXY-SND Infant foam sandal, for use with OxyTip+™ Sensitive Skin sensor 3/box

Monitors



S/5, AS/3, CS/3 Modular Monitors with M-series modules


Cardiocap/5, Cardiocap

S/5 Light Monitor

Capnomac Ultima

Satlite Plus II

Connecting Cables

 Interconnect cable with "N" (Datex) connector and "UN" connector (1/box)
OXY-SL3 (3m)

 Integrated cable with "N" connector and sensor

Sensor Type

OXY-F-UN Finger sensor with UN connector, 1 m, requires interconnect cable (1/box)

OXY-W-UN Wrap sensor with UN connector, 1 m, requires interconnect cable (1/box)

OXY-SE-3 Sensitive Skin sensor with UN connector, 4 m, requires interconnect cable (3/box)

OXY-AP-25 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (25/box)

OXY-AP-10 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (10/box)

OXY-AF-10 AllFit adhesive sensor, 0.9 m, extra adhesive tape included, for adult, pediatric and neonatal use, requires interconnect cable (10/box)

OXY-E4-N Ear sensor with N (Datex) connector, 4 m, headband and clip included (1/box)

OXY-F4-N Finger sensor with N (Datex) connector, 4 m, (1/box)

Accessories

- OXY-RTW** OxyTip+™ wide replacement tape, adhesive, for use with OxyTip+™ Wrap, OxyTip+™ Sensitive Skin sensor or OxyTip+™ AllFit 100/box
- OXY-RWL** Foam wrap replacement, large, weight range ≥3 kg, for use with OxyTip+™ Wrap sensor 24/box
- OXY-RWM** Foam Wrap replacement, medium, weight range ≥3 kg, for use with OxyTip+™ Sensitive Skin Sensor 24/box
- OXY-RWS** Foam Wrap replacement, small, weight range <3 kg or small finger, for use with OxyTip+™ Sensitive Skin sensor 24/box
- OXY-RTB** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, bears 100/box
- OXY-RT** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, blue 100/box
- OXY-SND** Infant foam sandal, for use with OxyTip+™ Sensitive Skin sensor 3/box

(See page 9 for more sensors)

Monitors



S/5, AS/3, CS/3 Modular Monitors with M-series modules


Cardicap/5, Cardicap

S/5 Light Monitor

Capnomac Ultima

Satlite Plus II

Connecting Cables

 Interconnect cable with "N" (Datex) connector and "DB" connector (1/box)
 OXY-C1 (1.5m)
 OXY-C3 (3m)
 OXY-C7 (7m)

 Integrated cable with "N" connector and sensor

Sensor Type

OXY-F-DB Finger sensor with DB connector, 2 m, requires interconnect cable with DB connector (1/box)

OXY-W-DB Wrap sensor with DB connector, 2 m, requires interconnect cable with "DB" connector (1/box)

OXY-E-DB Ear sensor with DB connector, 2 m, requires interconnect cable with "DB" connector (1/box)

OXY-E4-N Ear sensor with N (Datex) connector, 4 m, headband and clip included (1/box)

OXY-F4-N Finger sensor with N (Datex) connector, 4 m, 1/box)

Accessories

- OXY-RTW** OxyTip+™ wide replacement tape, adhesive, for use with OxyTip+™ Wrap, OxyTip+™ Sensitive Skin sensor or OxyTip+™ AllFit 100/box
- OXY-RWL** Foam wrap replacement, large, weight range ≥3 kg, for use with OxyTip+™ Wrap sensor 24/box
- OXY-RWM** Foam Wrap replacement, medium, weight range ≥3 kg, for use with OxyTip+™ Sensitive Skin Sensor 24/box
- OXY-RWS** Foam Wrap replacement, small, weight range <3 kg or small finger, for use with OxyTip+™ Sensitive Skin sensor 24/box
- OXY-RTB** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, bears 100/box
- OXY-RT** OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, blue 100/box
- OXY-SND** Infant foam sandal, for use with OxyTip+™ Sensitive Skin sensor 3/box

(See page 7 for more sensors)

Monitors



S/5, i4 Modular Monitors with E-series modules


Corometrics® 259 with TruSignal™
Ohmeda Technology

Connecting Cables



Interconnect cable, Ohmeda
OXY-ES3 (3m)

Sensor Type



OXY-F-UN Finger sensor with UN connector, 1 m, requires interconnect cable (1/box)

OXY-W-UN Wrap sensor with UN connector, 1 m, requires interconnect cable (1/box)

OXY-E-UN Ear sensor with UN connector, 1 m, requires interconnect cable (1/box)

OXY-SE-3 Sensitive Skin sensor with UN connector, 4 m, requires interconnect cable (3/box)

OXY-AP-25 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (25/box)

OXY-AP-10 Adult/Pediatric adhesive sensor, for adult and pediatric use, requires interconnect cable (10/box)

OXY-AF-10 AllFit adhesive sensor, 0.9 m, extra adhesive tape included, for adult, pediatric and neonatal use, requires interconnect cable (10/box)

Accessories

- OXY-RTW OxyTip+™ wide replacement tape, adhesive, for use with OxyTip+™ Wrap, OxyTip+™ Sensitive Skin sensor or OxyTip+™ AllFit 100/box
- OXY-RWL Foam wrap replacement, large, weight range ≥3 kg, for use with OxyTip+™ Wrap sensor 24/box
- OXY-RWM Foam Wrap replacement, medium, weight range ≥3 kg, for use with OxyTip+™ Sensitive Skin Sensor 24/box
- OXY-RWS Foam Wrap replacement, small, weight range <3 kg or small finger, for use with OxyTip+™ Sensitive Skin sensor 24/box
- OXY-RTB OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, bears 100/box
- OXY-RT OxyTip+™ replacement tape, for use with OxyTip+™ AllFit Sensor, blue 100/box
- OXY-SND Infant foam sandal, for use with OxyTip+™ Sensitive Skin sensor 3/box

Monitors



GE MRI Monitor



S/5 MRI Monitor

Integrated cable and sensor type



OXY-FMR

Integrated finger MRI sensor for adults, rectangle connector (6 m)



OXY-WMR

Integrated wrap MRI sensor for pediatrics, rectangle connector (6 m)



SAS-FMR

Integrated finger MRI sensor for adults, round connector (6 m)



SAS-WMR

Integrated wrap MRI sensor for pediatrics, round connector (6 m)

Monitors



Dash 3000/4000/5000 with Masimo option
 Solar 8000/8000M, i/9500 with Tram Module 451M/851M
 CardioLab®/Mac-Lab® with Tram Module 451M
 Corometrics® 259

Connecting Cables



Interconnect cable, Masimo LNOP, SpO₂
 2017002-001 (3.6m)
 2017002-003 (2.4m)

Sensor Type



- 2002799-001 Masimo LNOP Reusable Finger Sensor LNOP/DCI Pediatric
- 2002800-001 Masimo LNOP Reusable Finger Sensor Adult, LNOP/DCI
- 2010463-001 Masimo LNOP Reusable Multisite Sensor LNOP-YI
- 2027274-001 Masimo LNOP Reusable Tip-Clip Ear Sensor LNOP TC-I
- 2009745-001 Masimo LNOP Reusable Finger Sensor Adult DC-195
sensor connects directly to the interconnect cable

Accessories

- 2010466-001 Masimo Replacement Posey Wrap, LNOP-NeoPt-L, Neonatal (12/box.)
- 2010467-001 Masimo Tape Bag for LNOP-NEO (100/box)
- 2010468-001 Masimo Tape Cleanshield Multisite, LNOP-YI (100/box)
- 2010469-001 Masimo Disposable Standard Multisite Wrap, Adult/Ped/Neonatal Adhesive Attachment Wraps for use with LNOP-YI Multisite Reusable Sensor (100/box)
- 2010470-001 Masimo Tape Standard Petite Wrap, LNOP-YI (100/box)
- 2010471-001 Masimo Adhesive Tape for LNOP-YI (12/box.)

(See page 17 for more sensors)

Monitors




Dash 3000/4000/5000 with Masimo option

Solar 8000/8000M, i/9500 with Tram Module 451M/851M

CardioLab®/Mac-Lab® with Tram Module 451M

Corometrics® 259












Connecting Cables



Interconnect cable, Masimo LNOP, SpO₂
 2017002-001 (3.6m)
 2017002-003 (2.4m)

Sensor Type

Sensors connect directly to the interconnect cable

- 
2010458-001 Masimo LNOP Disposable Adhesive Sensor LNOP-ADT, Adult (20/box)
- 
2010459-001 Masimo LNOP Disposable Adhesive Sensor LNOP-PDT, Pediatric (20/box)
- 
2010461-001 Masimo LNOP Disposable Adhesive Sensor LNOP-NeoPt, Neonatal (20/box)
- 
2010460-001 Masimo LNOP Disposable Adhesive Sensor Bridge LNOP-NEO, (20/box)
- 
2017089-001 Masimo LNOP Disposable Adhesive Sensor LNOP-Neo-L, Neonatal (20/box)
- 
2017090-001 Masimo LNOP Disposable Adhesive Sensor LNOP-NeoPt-L, Neonatal (20/box)
- 
2027269-001 Masimo LNOP Disposable Adhesive Sensor Transparent Tape LNOP, Adult (20/box)
- 
2027270-001 Masimo LNOP Disposable Adhesive Sensor Transparent Tape LNOP, Pediatric (20/box)
- 
2027272-001 Masimo LNOP Disposable LNOP Hi Fi Sensor Neonatal/Adult (20/box)
- 
2027271-001 Masimo LNOP Disposable LNOP Hi Fi Sensor Infant/Pediatric (20/box)
- 
2027273-001 Masimo LNOP Disposable LNOP Blue Infant Thumb/Toe Sensor (10/box)

Accessories

- 2010466-001 Masimo Replacement Posey Wrap, LNOP-NeoPt-L, Neonatal (12/box.)
- 2010467-001 Masimo Tape Bag for LNOP-NEO (100/box)
- 2010468-001 Masimo Tape Cleanshield Multisite, LNOP-YI (100/box)
- 2010469-001 Masimo Disposable Standard Multisite Wrap, Adult/Ped/Neonatal Adhesive Attachment Wraps for use with LNOP-YI Multisite Reusable Sensor (100/box)
- 2010470-001 Masimo Tape Standard Petite Wrap, LNOP-YI (100/box)
- 2010471-001 Masimo Adhesive Tape for LNOP-YI (12/box.)

(See page 15 for more sensors)

Monitors



DINAMAP™ ProCare with Masimo option


DINAMAP™ ProSeries V2 with Masimo option

DINAMAP™ Pro 1000 V3 with Masimo option

Connecting Cables



Masimo PC08 LNOP interconnect cable
2009743-001 (2.4m)



Masimo PC04 LNOP interconnect cable
2027278-001 (1.2m)



Integrated Cable

Sensor Type



2002799-001 Masimo LNOP Reusable Finger Sensor LNOP/DCI Pediatric



2002800-001 Masimo LNOP Reusable Finger Sensor Adult, LNOP/DCI



2010463-001 Masimo LNOP Reusable Multisite Sensor LNOP-YI



2027274-001 Masimo LNOP Reusable Tip-Clip Ear Sensor LNOP TC-I



2009745-001 Masimo LNOP Reusable Finger Sensor Adult DC-195
sensor connects directly to the interconnect cable



2027276-001 Masimo LNOP DCSC Adult Reusable Sensor with 0.9m cable



2027277-001 Masimo LNOP DC-I-DC-12 Reusable Adult Sensor with 3.6m cable

Accessories

- 2010466-001 Masimo Replacement Posey Wrap, LNOP-NeoPt-L, Neonatal (12/box.)
- 2010467-001 Masimo Tape Bag for LNOP-NEO (100/box)
- 2010468-001 Masimo Tape Cleanshield Multisite, LNOP-YI (100/box)
- 2010469-001 Masimo Disposable Standard Multisite Wrap, Adult/Ped/Neonatal Adhesive Attachment Wraps for use with LNOP-YI Multisite Reusable Sensor (100/box)
- 2010470-001 Masimo Tape Standard Petite Wrap, LNOP-YI (100/box)
- 2010471-001 Masimo Adhesive Tape for LNOP-YI (12/box.)

(See page 21 for more sensors)

Monitors




DINAMAP™ ProCare with Masimo option

DINAMAP™ ProSeries V2 with Masimo option

DINAMAP™ Pro 1000 V3 with Masimo option

Connecting Cables



Masimo PC08 LNOP interconnect cable
2009743-001 (2.4m)












Masimo PC04 LNOP interconnect cable
2027278-001 (1.2m)

Accessories

- 2010466-001 Masimo Replacement Posey Wrap, LNOP-NeoPt-L, Neonatal (12/box.)
- 2010467-001 Masimo Tape Bag for LNOP-NEO (100/box)
- 2010468-001 Masimo Tape Cleanshield Multisite, LNOP-YI (100/box)
- 2010469-001 Masimo Disposable Standard Multisite Wrap, Adult/Ped/Neonatal Adhesive Attachment Wraps for use with LNOP-YI Multisite Reusable Sensor (100/box)
- 2010470-001 Masimo Tape Standard Petite Wrap, LNOP-YI (100/box)
- 2010471-001 Masimo Adhesive Tape for LNOP-YI (12/box.)

Sensor Type

Sensors connect directly to the interconnect cable

-  2010458-001 Masimo LNOP Disposable Adhesive Sensor LNOP-ADT, Adult (20/box)
-  2010459-001 Masimo LNOP Disposable Adhesive Sensor LNOP-PDT, Pediatric (20/box)
-  2010461-001 Masimo LNOP Disposable Adhesive Sensor LNOP-NeoPt, Neonatal (20/box)
-  2010460-001 Masimo LNOP Disposable Adhesive Sensor Bridge LNOP-NEO, (20/box)
-  2017089-001 Masimo LNOP Disposable Adhesive Sensor LNOP-Neo-L, Neonatal (20/box)
-  2017090-001 Masimo LNOP Disposable Adhesive Sensor LNOP-NeoPt-L, Neonatal (20/box)
-  2027269-001 Masimo LNOP Disposable Adhesive Sensor Transparent Tape LNOP, Adult (20/box)
-  2027270-001 Masimo LNOP Disposable Adhesive Sensor Transparent Tape LNOP, Pediatric (20/box)
-  2027272-001 Masimo LNOP Disposable LNOP Hi Fi Sensor Neonatal/Adult (20/box)
-  2027271-001 Masimo LNOP Disposable LNOP Hi Fi Sensor Infant/Pediatric (20/box)
-  2027273-001 Masimo LNOP Disposable LNOP Blue Infant Thumb/Toe Sensor (10/box)

(See page 19 for more sensors)

Monitors



Dash 2000/3000/4000 with Nellcor option

Dash 3000/4000/5000 with Ohmeda option

Solar 8000/8000M,i/9500 with Tram Module 451/451N/851/851N

CardioLab®/Mac-Lab® with Tram Module 451/451N

Corometrics® 120/2120 series

Eagle 1000/3000/4000

Responder 3000 with Nellcor option


Dash 3000/4000/5000 with Nellcor OxiMax option

Solar 8000/8000M, i9500 with Tram module 451N5/851N5

Corometrics® 250 Series with Nellcor OxiMax option



Connecting Cables



Interconnect cable Nellcor Oxismart
2006644-001 (2.9m)



Interconnect cable Nellcor OxiMax
2021406-001 (3m)



Interconnect cable Nellcor Oximax for Corometrics only
2025350-001 (3m)

Sensor Type



70124021 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)



70124033 Nellcor SpO₂ Multisite Sensor D-YS



70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)



70124035 Nellcor SpO₂ Sensor OXY-A/N, Neonatal/Adult (1/box)



414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)



70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)



70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)



70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)



70124032 Nellcor SpO₂ Sensor MAX-N, Adhesive, Neonatal (24/box)



407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)



2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)

Accessories

- 2016130-001 Nellcor Tape ADH-A/N, for use with 70124035 (100/box)
- 2016131-001 Nellcor Tape ADH-P/I for use with Oxi-P/I Sensors (100/box)

Monitors



DINAMAP™ ProCare with Nellcor OxiMax option

DINAMAP™ ProSeries V2 with Nellcor OxiMax option


DINAMAP™ Pro 1000 V3 with Nellcor OxiMax option

Connecting Cables



Interface cable Nellcor OxiMax
DOC-10
2008773-001 (2.9m)

Sensor Type



- 70124021 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)
- 70124033 Nellcor SpO₂ Multisite Sensor D-YS
- 70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)
- 70124035 Nellcor SpO₂ Sensor OXY-A/N, Neonatal/Adult (1/box)
- 414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)
- 70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)
- 70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)
- 70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)
- 70124032 Nellcor SpO₂ Sensor MAX-N, Adhesive, Neonatal (24/box)
- 407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)
- 2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)

Accessories

- 2016130-001 Nellcor Tape ADH-A/N, for use with 70124035 (100/box)
- 2016131-001 Nellcor Tape ADH-P/I for use with Oxi-P/I Sensors (100/box)

Monitors

DINAMAP™ ProSeries V1 and DINAMAP™ Pro 1000 V2 with Nellcor option



Connecting Cables



Interconnect cable Nellcor SCP-10
SCP-10 (2.9 m)



Interconnect cable Nellcor GE Corometrics® 118 Monitor
4033CAX (3 m)



Interconnect cable Nellcor EC8 (2.4 m)



Interconnect cable Nellcor 30344358 (2.9 m)



Interconnect cable Nellcor OxiMax 2026751-001 (2.9 m)

Sensor Type



70124021 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)



70124033 Nellcor SpO₂ Multisite Sensor D-YS



70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)



70124035 Nellcor SpO₂ Sensor OXI-A/N, Neonatal/Adult (1/box)



414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)



70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)



70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)



70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)



70124032 Nellcor SpO₂ Sensor MAX-N, Adhesive, Neonatal (24/box)



407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)



2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)

Accessories

- 2016130-001 Nellcor Tape ADH-A/N, for use with 70124035 (100/box)
- 2016131-001 Nellcor Tape ADH-P/I for use with Oxi-P/I Sensors (100/box)

Monitors



M-NSAT Modules for AS/3, CS/3, S/5 Modular Monitors

Connecting Cables




Interconnect cable Nellcor MC-10
903131 (3m)



Interconnect cable Nellcor OxiMax DOC-10
2008773-001 (2.9m)

Sensor Type



- 70124021 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)
- 70124033 Nellcor SpO₂ Multisite Sensor D-YS
- 70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)
- 70124035 Nellcor SpO₂ Sensor OXI-A/N, Neonatal/Adult (1/box)
- 414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)
- 70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)
- 70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)
- 70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)
- 70124032 Nellcor SpO₂ Sensor MAX-N, Adhesive, Neonatal (24/box)
- 407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)
- 2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)



E-NSAT Modules for i4 Modular Monitors

Accessories

- 2016130-001 Nellcor Tape ADH-A/N, for use with 70124035 (100/box)
- 2016131-001 Nellcor Tape ADH-P/I for use with Oxi-P/I Sensors (100/box)

Monitors



Dash 2000/3000/4000 with Nellcor option
 Dash 3000/4000/5000 with Ohmeda option
 Solar 8000/8000M, i/9500 with Tram Module 451/451N/851/851N
 Eagle 1000 (rectangle connector)
 Responder 3000 with Nellcor option
 Corometrics 120/2120 series

DINAMAP™ ProSeries V1 and DINAMAP™ Pro 1000 V2 with Nellcor option



Cardioserv Eagle 1000 (round connector)

Corometrics 118/556

Connecting Cables



Interconnect cable Nellcor Oxismart 2006644-001 (2.9m)



Interconnect cable Nellcor SCP-10 SCP-10 (2.9 m)



Interconnect cable Nellcor 30344358 (2.9 m)



Interconnect cable Nellcor GE Corometrics® 118 Monitor 4033CAX (3 m)


Sensor Type



2023211-001 Nellcor Compatible Finger Sensor with DB9 connector 1M



2023215-001 Nellcor Compatible Ear Sensor with DB9 connector 1M



2023216-001 Nellcor Compatible Neonatal-Adult Adhesive Sensor with DB9 connector (24/box)



2023217-001 Nellcor Compatible Adult/Pediatric Adhesive Sensor DB9 connector (24/box)

Monitors



DINAMAP™ ProCare with Masimo option
 DINAMAP™ ProSeries V2 with Masimo option
 DINAMAP™ Pro 1000 V3 with Masimo option



Dash 3000/4000/5000 with Masimo option
 Solar 8000/8000M, i/9500 with Tram Module 451M/851M
 CardioLab®/Mac-Lab® with Tram Module 451M
 Corometrics® 259

Connecting Cables



Masimo LNC-10 interconnect cable for GE Dinamap
 2027263-001 (3m)



Masimo LNC-10 interconnect cable for Dash, Tram, Corometrics
 2027263-002 (3m)

Sensor Type

Reusable Sensors

- 2027258-001 Masimo LNCS DC-I Reusable Adult Sensor
- 2027259-001 Masimo LNCS DC-IP Reusable Pediatric Sensor
- 2027261-001 Masimo LNCS TC-I TipClip Reusable Ear Sensor

Adhesive Sensors

- 2027253-001 Masimo LNCS Adult Adhesive Sensor (20/box)
- 2027254-001 Masimo LNCS Pediatric Adhesive Sensor (20/box)
- 2027255-001 Masimo LNCS Infant Adhesive Sensor (20/box)
- 2027256-001 Masimo LNCS Neonatal Adhesive Sensor (20/box)
- 2027257-001 Masimo LNCS Neonatal PT Adhesive Sensor (20/box)